


# Materials and Methods

## RNA World and the Origin of life

### 1. RNA Production



General scheme of RNA production for P1-Loop, T<sub>2</sub>S1, substrate and ribozymes sequences

*In vitro* replication of main sequences of sunY and Azoarcus mutants were carried out in an Eppendorf Tube, using the following protocol in quantities and the later the PCR program used.

All sequences were preserved at -20°C before use.

	Stock solution	Unit	V(μl)	Final conc.	Unit
Ribozyme template	1	ng/μL	10	0,2	ng/μL
Buffer	5	X	10	1	X
Forward primer	10	μM	3,75	0,75	μM
Reverse primer	10	μM	3,75	0,75	μM
dNTP Mix	10	mM	1	0,2	mM
Polymerase (Phusion)	2	U/μL	0,5	0,02	U/μL
Water	-	-	21	50	
		<b>Total</b>	<b>50</b>		

PCR program:			
Step 1	98 °C	30 sec	Denaturation
Step 2	98 °C	10 sec	
Step 3	59 °C	30 sec	Annealing
Step 4	72 °C	30 sec	Elongation
Step 5	Go to step 2	30 times	
Step 6	72 °C	3 min	

After PCR, 2μL of each PCR product + 1μl stopping buffer were used for the agarose gel analysis. 25 mL of 2% gel were prepared by mixing 1 tablet of TopVision Agarose Tablets, containing 0,5 g agarose, and 25 mL TAE buffer 1X. Solution was heated in a microwave until it appeared clean. The gel was pre-run 10 minutes using the same buffer, at 12 V/cm gel was left 25 min. The purpose of analysis was to determine if the PCR cycles had been successful in replicating the sequence.

Only then RNA production can begin with the usage of HiScribe™ T7 High Yield RNA Synthesis Kit which uses the T7 RNA Polymerase. Protocol was followed as indicated by the fabricant each reaction was carried out in 20 μL. After the reaction purification was needed. It comprised several steps:

#### Phenol-chloroform extraction in acid medium:

Resuspension of the IVT solution with water (nuclease free 100μl) and transfer the solution into a 1,5 mL Eppendorf tube, under the hood 100 μL of phenol-chloroform pH 4 were add, vortex a few minutes, then centrifuge 4 minutes at 11 000 rpm, during the centrifugation, add 0.1 volume of 3M NaAc (sodium acetate) in new tubes at the end of the centrifugation, gently pipette the upper aqueous phase (containing RNAs) into the tubes containing NaAc add 2.5 volumes of cold 100% ethanol (stored at -20°C) and mix by up and down (the solution should become cloudy, RNAs precipitate quickly) let the RNAs precipitate at -20°C for 1h or overnight, centrifuge for 1h at 4°C 13 000 rpm, remove the supernatant, a white pellet should be visible. Wash the pellet with 150 μL of cold 70% ethanol (stored at -20°C) very gently (without dissolving the pellet) then centrifuge briefly and remove the supernatant. Repeat the wash a second time removing as much ethanol as possible cover the tube with parafilm (to prevent the pellet from escaping) and pierce some holes with a syringe dry the pellet in the concentrator under vacuum for a few minutes resuspend the pellet with 40 μL of water (pipette well until the pellet is dissolved)

## DNase I Treatment

We begin with 5  $\mu\text{L}$  of DNase I buffer (10X) and add 5  $\mu\text{L}$  of enzyme and incubate 15 min at 37°C immediately add 50  $\mu\text{L}$  of stopping buffer. The sample can be stored at -20°C or can be deposited directly on extraction gel RNA extraction on acrylamide (PAGE) gel. Run the 8% acrylamide (PAGE) gel and preheat the gel (24W, 15 minutes). Denature the samples by heat treatment in thermocycler (5 minutes at 75°C) migration is done until the bottom marker is at the bottom of the gel (about 3 hours, 24W). Remove the gel, cut the top right edge to make a marker, and place it between 2 plastic films, place the gel on a silicon plate, illuminate at 254 nm and mark the outline with a marker, cut the gel with a razor blade. Place the gel pieces in the sample tubes, crush with a 1mL cone the gel on the wall of the tube and add 0.5 mL of 0.3M NaAc. Incubation in Thermomix for 5 hours at 26°C, finally centrifugation at 450-500 rpm.

## RNA purification by EtOH precipitation

Transfer the contents of the tubes (taking care not to over-pipette the gel ends) into COSTAR columns and centrifuge for 4 minutes at 11 000 rpm (-add 1 $\mu\text{L}$  of glycogen if necessary). Add 2.5 volumes of 100% EtOH cooled and incubate at -20°C for 1 hour. Centrifuge at 13 000 rpm at 4°C for 1 hour. Wash the pellets with 150 $\mu\text{L}$  of cold 70% EtOH. Evaporate the ethanol (under vacuum or hood), resuspend in 20 $\mu\text{L}$  of water and measure the concentration with the nanodrop.

## 2. Catalytic Assays

All catalytic assays were performed in the same fashion manner. According to the variability of experimental conditions several catalytic assays could be run with the same controls. All sequences were introduced leaving at the end the ribozyme when the start point was taken immediately. All Eppendorf tubes were kept in a thermostat at 40 °C for minutes. Only for the actual assay samples (2  $\mu\text{L}$ ) were taken at 0, 2, 4, 7 and 24 hours for SunY and 0, 0,5, 1, 2, 4, 7 and 24 hours for Azoarcus. For controls only initial and final points were acquired. These samples were left in 6  $\mu\text{L}$  stopping buffer (10 mM EDTA, 80% formamide, 5% bromophenol blue, 5% xylene cyanol). Before loading of the samples to the gel denaturation by temperature was done in thermocycler at 80°C for 5 minutes.

Samples were analyzed in 18% PAGE (Polyacrylamide gel with 50% as denaturing agent). 40 mL of Polyacrylamide solution were prepared using 198  $\mu\text{L}$  10% APC and 66  $\mu\text{L}$  TEMED as polymerization agents. 4  $\mu\text{L}$  of sample were introduced by lane. 4  $\mu\text{L}$  ladders of 1, 0,5 and 0,25 mM were used of molecular weights of 21, 25, 36, 46 and 55 bp for all cases except in P1-Loop-Azo ligation where 25 bp standard was replace by a 26 bp standard. Gel was pre-run for 10 min in TBE buffer 1X, at 25 V/cm gel was left around an hour. Visualization of the gel was done using SYBR Gold system with a Blue Tray.

### Protocol for P1 Loop ligation (Azoarcus vs SunY experiment)

Catalytic assay				
Reagent	Stock sol.	Units	V( $\mu\text{L}$ )	Final conc.
Ribozyme	16	$\mu\text{M}$	0,3	0,5
P1 Loop	100	$\mu\text{M}$	1,5	15
S1_LL	100,0	$\mu\text{M}$	1,5	15
MgCl2	240	mM	2,5	60
Water	up to		4,2	
		Sum	10	

Ribozyme + P1 Loop Control				
Reagent	Stock sol.	Units	V( $\mu\text{L}$ )	Final conc.
Ribozyme	16	$\mu\text{M}$	0,3	0,5
P1 Loop	100	$\mu\text{M}$	1,2	15
S1_LL	100,0	$\mu\text{M}$	0,0	0
MgCl2	240	mM	2,0	60
Water	up to		4,6	
		Sum	8	

P1 Loop + S1 Control				
Reagent	Stock sol.	Units	V( $\mu\text{L}$ )	Final conc.
Ribozyme	16	$\mu\text{M}$	0,0	0
P1 Loop	100	$\mu\text{M}$	1,2	15
S1_LL	100,0	$\mu\text{M}$	1,2	15
MgCl2	240	mM	2,0	60
Water	up to		3,6	
		Sum	8	

Ribozyme + S1 Control				
Reagent	Stock sol.	Units	V( $\mu\text{L}$ )	Final conc.
Ribozyme	16	$\mu\text{M}$	0,3	0,5
P1 Loop	100	$\mu\text{M}$	0,0	0
S1_LL	100,0	$\mu\text{M}$	1,2	15
MgCl2	240	mM	2,0	60
Water	up to		4,6	
		Sum	8,0	

P1 Loop Control				
Reagent	Stock sol.	Units	V( $\mu\text{L}$ )	Final conc.
Ribozyme	16	$\mu\text{M}$	0,0	0
P1 Loop	100	$\mu\text{M}$	1,2	15
S1_LL	100,0	$\mu\text{M}$	0,0	0
MgCl2	240	mM	2,0	60
Water	up to		4,8	
		Sum	8	

S1 Control				
Reagent	Stock sol.	Units	V( $\mu\text{L}$ )	Final conc.
Ribozyme	16	$\mu\text{M}$	0,0	0
P1 Loop	100	$\mu\text{M}$	0,0	0
S1_LL	100,0	$\mu\text{M}$	1,2	15
MgCl2	240	mM	2,0	60
Water	up to		4,8	
		Sum	8	

Ribozyme Control				
Reagent	Stock sol.	Units	V( $\mu\text{L}$ )	Final conc.
Ribozyme	16	$\mu\text{M}$	0,3	0,5
P1 Loop	100	$\mu\text{M}$	0,0	0
S1_LL	100,0	$\mu\text{M}$	0,0	0
MgCl2	240	mM	2,0	60
Water	up to		5,8	
		Sum	8	

P1-Loop and T\_2S1 analysis with A and B *Azoarcus* modifications were done in similar manner 10  $\mu$ L for catalytic assays and 8  $\mu$ L for controls. Only concentrations were modified. Assay **A-B-P1-Loop**: Substrates 15  $\mu$ M. Assay **A-B-P1-Loop**: Substrate S1 15  $\mu$ M + 2.5  $\mu$ M T\_2S1.

### Protocol for SunY Mg optimization in Water

MgCl2 50 mM				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,3	0,5
T_2S_RNA	100	$\mu$ M	0,5	5
S1_LL	100,0	$\mu$ M	1,5	15
MgCl2	240	mM	2,1	50
Water	up to		5,6	
Sum			10	

MgCl2 75 mM				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,3	0,5
T_2S_RNA	100	$\mu$ M	0,5	5
S1_LL	100,0	$\mu$ M	1,5	15
MgCl2	240	mM	3,1	75
Water	up to		4,6	
Sum			10,0	

MgCl2 100 mM				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,3	0,5
T_2S_RNA	100	$\mu$ M	0,5	5
S1_LL	100,0	$\mu$ M	1,5	15
MgCl2	240	mM	4,2	100
Water	up to		3,5	
Sum			10	

Sun Y + T_2S Control				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,3	0,5
T_2S_RNA	100	$\mu$ M	0,4	5
S1_LL	100,0	$\mu$ M	0,0	0
MgCl2	240	mM	3,3	100
Water	up to		4,0	
Sum			8	

T_2S + S1 Control				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,0	0
T_2S_RNA	100	$\mu$ M	0,4	5
S1_LL	100,0	$\mu$ M	1,2	15
MgCl2	240	mM	3,3	100
Water	up to		3,1	
Sum			8	

Sun Y + S1 Control				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,3	0,5
T_2S_RNA	100	$\mu$ M	0,0	0
S1_LL	100,0	$\mu$ M	1,2	15
MgCl2	240	mM	3,3	100
Water	up to		3,2	
Sum			8,0	

T_2S Control				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,0	0
T_2S_RNA	100	$\mu$ M	0,4	5
S1_LL	100,0	$\mu$ M	0,0	0
MgCl2	240	mM	3,3	100
Water	up to		4,3	
Sum			8	

S1 Control				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,0	0
T_2S_RNA	100	$\mu$ M	0,0	0
S1_LL	100,0	$\mu$ M	1,2	15
MgCl2	240	mM	3,3	100
Water	up to		3,5	
Sum			8	

SunY Control				
Reagent	Stock sol.	Units	V( $\mu$ l)	Final conc.
SunY	16	$\mu$ M	0,3	0,5
T_2S_RNA	100	$\mu$ M	0,0	0
S1_LL	100,0	$\mu$ M	0,0	0
MgCl2	240	mM	3,3	100
Water	up to		4,4	
Sum			8	

### Protocol for SunY Mg optimization in EPPS 30 mM

Volume 10 $\mu$ L			MgCl2 0 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	0,0	0	mM
Buffer	200	mM	1,5	30	mM
Water	up to		5,9		
Sum			10,0		

Volume 10 $\mu$ L			MgCl2 25 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	1,0	25	mM
Buffer	200	mM	1,5	30	mM
Water	up to		4,9		
Sum			10,0		

Volume 10 $\mu$ L			MgCl2 50 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	2,1	50	mM
Buffer	200	mM	1,5	30	mM
Water	up to		3,9		
Sum			10,0		

Volume 10 $\mu$ L			MgCl2 75 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	3,1	75	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,8		
Sum			10,0		

Volume 10 $\mu$ L			MgCl2 100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	4,2	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		1,8		
Sum			10,0		

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,4	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,4	5	$\mu$ M
S1_LL	100,0	$\mu$ M	0,0	0	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,3		
Sum			8,0	Sun Y + P1	

CONTROLS

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,4	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,0	0	$\mu$ M
S1_LL	100,0	$\mu$ M	1,2	15	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		1,5		
Sum			8,0	Sun Y + S1	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,0	0	$\mu$ M
T-2S	100	$\mu$ M	0,4	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,2	15	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		1,6		
Sum			8,0	S1 + P1	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,4	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,0	0	$\mu$ M
S1_LL	100,0	$\mu$ M	0,0	0	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,7		
Sum			8,0	Sun Y	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,0	0	$\mu$ M
T-2S	100	$\mu$ M	0,4	5	$\mu$ M
S1_LL	100,0	$\mu$ M	0,0	0	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,8		
Sum			8,0	P1	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,0	0	$\mu$ M
T-2S	100	$\mu$ M	0,0	0	$\mu$ M
S1_LL	100,0	$\mu$ M	1,2	15	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,0		
Sum			8,0	S1	

### Protocol for SunY Mg optimization in Tris-HCl 30 mM

Volume			MgCl2		
10 $\mu$ L			0 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	0,0	0	mM
Buffer	200	mM	1,5	30	mM
Water	up to		5,9		
Sum			10,0		

Volume			MgCl2		
10 $\mu$ L			25 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	1,0	25	mM
Buffer	200	mM	1,5	30	mM
Water	up to		4,9		
Sum			10,0		

Volume			MgCl2		
10 $\mu$ L			50 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	2,1	50	mM
Buffer	200	mM	1,5	30	mM
Water	up to		3,9		
Sum			10,0		

Volume			MgCl2		
10 $\mu$ L			75 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	3,1	75	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,8		
Sum			10,0		

Volume			MgCl2		
10 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,6	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,5	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,5	15	$\mu$ M
MgCl2	240	mM	4,2	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		1,8		
Sum			10,0		

CONTROLS

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,4	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,4	5	$\mu$ M
S1_LL	100,0	$\mu$ M	0,0	0	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,3		
Sum			8,0	Sun Y + P1	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,4	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,0	0	$\mu$ M
S1_LL	100,0	$\mu$ M	1,2	15	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		1,5		
Sum			8,0	Sun Y + S1	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,0	0	$\mu$ M
T-2S	100	$\mu$ M	0,4	5	$\mu$ M
S1_LL	100,0	$\mu$ M	1,2	15	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		1,6		
Sum			8,0	S1 + P1	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,4	0,5	$\mu$ M
T-2S	100	$\mu$ M	0,0	0	$\mu$ M
S1_LL	100,0	$\mu$ M	0,0	0	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,7		
Sum			8,0	Sun Y	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,0	0	$\mu$ M
T-2S	100	$\mu$ M	0,4	5	$\mu$ M
S1_LL	100,0	$\mu$ M	0,0	0	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,8		
Sum			8,0	P1	

Volume			MgCl2		
8 $\mu$ L			100 mM		
name	stock	unit	vol( $\mu$ l)	final conc.	unit
SunY	9	$\mu$ M	0,0	0	$\mu$ M
T-2S	100	$\mu$ M	0,0	0	$\mu$ M
S1_LL	100,0	$\mu$ M	1,2	15	$\mu$ M
MgCl2	240	mM	3,3	100	mM
Buffer	200	mM	1,5	30	mM
Water	up to		2,0		
Sum			8,0	S1	

### 3. Image Processing

Processing was done in GelAnalyzer. Lanes were selected and use of the Ladder was used to determine the molecular weight of different fragments. Background estimation was done with rolling ball method and 5% pic fit. Deviation for error bars were done considering 10% of area outside of gaussian pic.